



- 580 Watts Cooling Capacity
- 7 L/min Pump Flow Rate
- Quiet Operation
- Accessible Cold Reservoir Tank
- -10°C to +30°C Operating Range
- Extremely Reliable

Circulate Externally

For applications requiring external re-circulation of cold water, the RC4 delivers a pressurised refrigerated water supply through either an open or closed loop system. The RC4 is ideal for cooling condenser systems, water jacketed devices, electrophoresis baths and heat exchangers, including those found in Ratek OM15C and OM25 incubator shakers. It's quiet operation makes it ideal for use in the laboratory, but it is equally suited to industrial cooling applications.



SPECIFICATIONS

| RC4 | |
|----------------------|---|
| Cooling capacity | Average 580 Watts between 30°C and 5°C |
| Reservoir dimensions | L275 x W200 x D200mm |
| Filling capacity | 12 litres including heat exchanger |
| Refrigerant type | R134a |
| Pump flow rate | 7 litres per minute at zero head, max pressure 3.2 psi |
| Condenser | ¼ hp air cooled hermetic condensing unit |
| Operating range | -10°C to +30°C (using suitable additives) |
| Thermostat | Adjustable low-temp safety cut-off (factory default is +4°C). For temperature control, the TH8000 is required. |
| Power input | 240VAC/690W |
| Overall dimensions | L480 x W340 x D510mm |
| Nett weight | 27 kg |

When used in a closed loop environment such as that found in a condenser system, the return water feed comes back directly under the pump's pressure. It is also possible to use the RC4 in an open-loop system such as an open water bath or container, where the return feed is not directly connected to the outlet. In this case the pump's suction then pulls the water back around the circuit making the RC4 ideal for refrigerating existing water-bath systems.

The RC4 includes an accessible 12 litre reservoir with lid that can accommodate a Ratek immersion heater circulator to provide a precision temperature controlled water supply which can be maintained between -10°C and +30°C (4°C to 30°C factory default). The TH8000 is ideal for this application.

The accessible reservoir can also be used directly as a small refrigerated bath with vessels placed directly into the tank, making the RC4 a very versatile and cost effective unit that can be used for a variety of applications around the lab.

Built to last with both ultra-reliable and quiet performance, the RC4 is made and supported in Australia and backed by Ratek's unparalleled warranty and after sales service.

RELATED PRODUCTS

| Product Code | Description |
|--------------|--|
| TH8000 | Precision immersion heater for accurate temperature controlled water supply. |

YOU MAY ALSO LIKE

RC1 - Immersion Cooler

The RC1 is a re-locatable immersion cooling coil that delivers up to 580 Watts of fixed cooling to virtually any bath or vessel.



AUSTRALIAN
MADE & SUPPORTED

3 YEAR
PARTS WARRANTY



Quality
ISO 9001
SAI GLOBAL