



- Repeatable Dial Speed Control
- Simple Operation
- Adjustable Tilt Angle
- 0-60 RPM Speed Adjustable
- 300 x 180mm Platform
- Can Create Multiple Mixing Actions

### Extremely Versatile

The GM1 is the perfect personal 3D rocker suitable for a wide range of mixing applications. The unique gyratory movement sometimes referred to as the "Belly Dancer" action tilts the platform in the horizontal plane through 360° to provide a gentle swirling motion, ideal for gentle mixing in a variety of different types of vessels. The compact design is ideal for mixing small quantities of vessels directly at your workstation.



### SPECIFICATIONS

GM1	
Mixing action	3D Gyratory (adjustable)
Platform size	300 x 180 mm
Speed range	0-60 RPM
Speed control type	Analogue
Motor	Swiss made DC motor
Transmission type	Gearbox
Power input	240VAC/20W
Overall dimensions	W300 x D190 x H220 mm
Nett weight	3 kg

The mixing action can quickly be adjusted by tilting the platform to the desired angle. The ability to adjust the tilt angle combined with a smooth precision speed control allows you to achieve the perfect mixing action for your application.

Ideal for DNA extractions and staining/de-staining procedures, plus low foaming agitation in blood and centrifuge tubes.

The GM1 is simple to operate yet offers precise repeatable speed control in a package that is built to last.

Backed by Ratek's unparalleled service, support and warranty, the GM1 is an essential tool for every workstation.

### RELATED PRODUCTS

Product Code	Description
AT1	Adhesive mat tray to suit GM1

### YOU MAY ALSO LIKE

#### ERPM4 - Medium Rocker

The ultra-reliable ERPM4 has a 350 x 350mm platform and simple yet repeatable dial speed control.



AUSTRALIAN  
MADE & SUPPORTED

3 YEAR  
PARTS WARRANTY



Quality  
ISO 9001  
SAI GLOBAL