



- Repeatably Digital Tacho Control
- 3 User Definable Program Presets
- Near Silent Operation
- 20-200 RPM Speed Adjustable
- 350 x 350mm Platform
- Large Range of Accessories Available

Precision Control

The RM2 couples a reciprocating (also known as linear or side to side) action with precision speed control and digital tachometer. It is ideal for creating a gentle washing action in plates, bags and flat-bottomed sample containers, but is equally suited to vigorous end to end shaking at higher speeds. The RM2 is perfectly suited for mixing Mojonnier flasks and the dedicated MR5 Mojonnier rack is available for this application.



SPECIFICATIONS

RM2	
Mixing action	Reciprocating (linear, side to side)
Platform size	350 x 350 mm
Speed range	20-200 RPM
Speed control type	Digital + Tacho
Timer	0-99 mins/hrs
LED speed & time display	Yes
User program presets	3
Motor	Long life super-quiet brushless AC
Transmission type	Silent Belt
Safety features	<ul style="list-style-type: none"> • Motion interrupt stop • Slow start • Unbalanced load stop
Power input	240V AC/100W
Overall dimensions	W370 x D400 x H160 mm
Nett weight	10.2 kg

With a super-quiet drive system incorporating slow acceleration to reduce spillage, the RM2 delivers smooth, precise and repeatable speed control via its large user-friendly digital control panel. The drive system also features out of balance and platform jam safety cut-out systems.

The brilliant red LED tachometer, timer and user program display ensures you know exactly what your samples are doing, visible from anywhere in the lab, whilst the audible program buzzer will alert you once a program is complete.

A wide range of accessories is available including universal racks, tulip clip trays, adhesive mat trays and dual-tier trays plus customised trays can be made on request.

RELATED PRODUCTS

Product Code	Description
AT5	Adhesive mat tray to suit medium mixers
RR5	Universal rack to suit medium mixers
TCT5	Tulip clip tray to suit medium mixers (clips sold separately)
DTT5	Dual-tier tray to suit medium mixers
MR5	Mojonnier rack tray to suit medium mixers

AUSTRALIAN
MADE & SUPPORTED

3 YEAR
PARTS WARRANTY



Quality
ISO 9001
SAI GLOBAL