



- Repeatabe Digital Tacho Control
- 3 User Definable Program Presets
- Adjustable Tilt Angle
- 4-57 RPM Speed Adjustable
- 350 x 350mm Platform
- Large Range of Accessories Available

Precision Control

The RPM5 is Ratek's most advanced rocking platform mixer and incorporates full precision digital control with our silent-belt drive system and a host of safety features. The platform tilt angle can also be adjusted to tailor the mixing action to suit your application from a gentle washing action right through to vigorous end to end mixing. The RPM5 is ideal for use with the MR5 Mojonier rack to provide gentle full end to end mixing of Mojonier flasks.



SPECIFICATIONS

RPM5	
Mixing action	Rocking with adjustable tilt angle
Platform size	350 x 350 mm
Speed range	4-57 RPM
Speed control type	Digital + Tacho
Timer	0-99 mins/hrs
LED speed & time display	Yes
User program presets	3
Motor	Long life super-quiet brushless AC
Transmission type	Silent Belt
Safety features	<ul style="list-style-type: none"> • Motion interrupt stop • Slow start • Unbalanced load stop
Power input	240VAC/100W
Overall dimensions	W350 x D400 x H200 mm
Nett weight	10.1 kg

With a super-quiet drive system incorporating slow acceleration to reduce spillage, the RPM5 delivers smooth, precise and repeatable speed control via it's large user-friendly digital control panel. The brilliant red LED display is visible from anywhere in the lab allowing the tachometer, timer and user programs to be monitored easily. The drive system also features out of balance and platform jam safety cut-out systems.

A wide range of accessories is available including universal racks, tulip clip trays and adhesive mat trays plus customised trays can be made on request.

When combined with the DTT5 dual-tier tray, the capacity of the mixer can be virtually doubled without affecting the benchtop footprint.

RELATED PRODUCTS

Product Code	Description
AT5	Adhesive mat tray to suit medium mixers
RR5	Universal rack to suit medium mixers
TCT5	Tulip clip tray to suit medium mixers (clips sold separately)
DTT5	Dual-tier tray to suit medium mixers
MR5	Mojonnier rack tray to suit medium mixers