

CAN BE
CUSTOM
MADE

RM4 Rotary Mixer - Digital



- Extremely Versatile
- Custom Made Carriers Available
- 6-57 RPM Speed Range
- End to End & Rolling Mixing Actions
- Digital Control & Large LED Displays
- Super Quiet Operation

The All-Round Rotary Mixer

The RM4 Rotary Mixer is very much an all-rounder, with near endless potential for different applications. Perfect for mixing larger laboratory bottles and flasks right through to wine and food products, the RM4 has made a home in a wide range of industries. Available with or without an inter-locked safety enclosure, the RM4 is easily customised to suit almost any job you can throw at it.



SPECIFICATIONS

RM4	
Mixing action	<ul style="list-style-type: none"> • End over end • Rolling (custom)
Capacity	4kgs evenly loaded
Speed range	6-57 RPM
Speed control type	Precision Digital + Tacho
Timer	0-99 mins or 0-99 hrs
LED speed and time display	Yes
Program presets	3 x user configurable presets
Motor and drive type	Long life super-quiet brushless AC with replaceable belt drive
Safety features	Motion interrupt cut out
Power input	240 VAC/130W
Overall dimensions	W455 x D400 x H335mm <ul style="list-style-type: none"> • Can be lengthened as required
Nett weight	10 kg

Supplied as standard, the RM4 can accommodate up to 4 x Ratek UMC1/2/3/4 clamps, allowing test tubes, vials, lab bottles, wine bottles, beer bottles or virtually any piece of laboratory glassware to be rotated to give an end to end mixing action. By mounting the clamps at 90°, a rolling action can be achieved inside the vessel on the same machine.

The RM4 can be supplied with custom designed carriers to suit your application, including fully enclosed hinged containers, and extended mixing bars with large numbers of clips to suit small vials or tubes.

The RM4 features precision digital speed control which is complemented by the Ratek silent-belt drive system and long life brushless motor, as well as large highly visible digital LED timer and tachometer displays and audible alarm.

Built in Melbourne, supported locally and supplied with our unparalleled warranty and backup service, the RM4 is the perfect fit for a wide range of applications, and makes light work of specialised mixing jobs.

RELATED PRODUCTS

Product Code	Description
SP707	Inter-locked safety enclosure for RM4
UMC1/2/3/4	Universal clamps to suit RM4

AUSTRALIAN
MADE & SUPPORTED

3 YEAR
PARTS WARRANTY



Quality
ISO 9001
SAI GLOBAL