

ASSEMBLY INSTRUCTIONS

For transportation the Retort stand and top tray are removed from the MTV1.

- Unpack the unit
- Fit the retort stand to the retaining block located on the rear panel with the flat to the back(the securing screws screw to the flat).
- Using the allen key provided tighten the socket screws to hold the retort stand in place.
- Slide the top plate onto the retort stand and use the clamp nut to hold the plate in the desired position.
- 4 rubber securing cups are supplied for additional stability, they are placed under the feet of the MTV1

ASSEMBLY INSTRUCTIONS

For transportation the Retort stand and top tray are removed from the MTV1.

- Unpack the unit
- Fit the retort stand to the retaining block located on the rear panel with the flat to the back(the securing screws screw to the flat).
- Using the allen key provided tighten the socket screws to hold the retort stand in place.
- Slide the top plate onto the retort stand and use the clamp nut to hold the plate in the desired position.
- 4 rubber securing cups are supplied for additional stability, they are placed under the feet of the MTV1

PRODUCT WARRANTY

Ratek equipment is covered by a 12 month parts and labour warranty effective from the date of purchase.

This warranty covers the repair or replacement of any parts or components found to be defective, subject to the service options listed below.

This warranty excludes any defect resulting from misuse, neglect, accidental damage, improper voltage, or any alteration which affects the performance of the equipment.

- Does not extend to any costs associated with the delivery, damage, or loss incurred during transport.
- Is in addition to any Statutory regulations and provisions implied by the Trade Practices Act and any relevant State or Federal Government obligations.

For warranty service the following options are available:

- Contact the supplier from where the equipment was purchased to arrange for service.
- Return the equipment direct to Ratek Instruments, 60 Wadhurst Drive, Boronia, Vic 3155.

Warranty repairs not returned to Ratek can be undertaken by authorised service personnel but are limited to a parts only warranty.

Note: Proof of purchase is required for all warranty repairs.



OPERATING INSTRUCTIONS

FOR

MTV1 Multi tube vortex mixer



FOR ADDITIONAL INFORMATION / TECHNICAL SUPPORT

TEL: 61 3 98872161

FAX: 61 3 98872163

EMAIL: sales@ratek.com.au